

Mouse Cell Cycle Experiments.

✧ Experimental preparation

1. BD Cycletest™ Plus DNA Kit
 - i. Buffer Solution (room temperature).
 - ii. Solution A (room temperature).
 - iii. Solution B (room temperature).
 - iv. Solution C (on ice, without light, PI).
2. Cell culture preparation
Cell amount: 3×10^5 cells/sample.

✧ Experimental procedure

1. Place cells from culture plate into tube.
2. Centrifuge the cell suspensions at 1500 rpm for 5 minutes at room temperature (20°C–25°C).
3. Aspirate the supernatant, leaving approximately 50 μ L of residual fluid in the tube to avoid disturbing the pellet.
4. Add 1 mL of Buffer Solution and resuspend the cells by gently vortexing at low speed.
5. Centrifuge for 5 minutes at 1500 rpm at room temperature (20°C–25°C).
6. Repeat step 3 through step 5.
7. Aspirate the supernatant, leaving approximately 50 μ L of residual fluid in the tube to avoid disturbing the pellet.
8. Resuspend the pellet in 1 mL of Buffer Solution by gently vortexing at low speed.
9. Count the cells.
10. Adjust the concentration to 1.0×10^6 cells/mL with Buffer Solution. (3×10^5 cells in the 300 μ L Buffer Solution)
This should provide sufficient cells for both a test sample tube and a control tube.
11. Centrifuge the cell suspensions at 1500 rpm for 5 minutes at room temperature (20°C–25°C).
12. Carefully decant all the supernatant.
13. Add 250 μ L of Solution A (trypsin buffer) to each tube and gently mix by tapping the tube by hand. ***Do not vortex.***
14. Incubate for 10 minutes at room temperature (20°C– 25°C). ***Do not aspirate Solution A.***
15. Add 200 μ L of Solution B (trypsin inhibitor and RNase buffer) to each tube and gently mix by tapping the tube by hand. ***Do not vortex.***
16. Incubate for 10 minutes at room temperature (20°C– 25°C). ***Do not aspirate Solution A and B.***
17. Add 200 μ L of cold (2°C–8°C) Solution C (PI stain solution) to each tube and gently mix by tapping the tube by hand. ***Do not vortex.***
18. Incubate for 10 minutes in the dark on ice or in the refrigerator (2°C–8°C).
19. Filter the sample into a tube.

The samples are now ready to be analyzed on the flow cytometer. Cap or cover the prepared tubes and store at 2°C– 8°C in the dark until flow cytometric analysis.