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Nature Protocols 17, 2431–2468 (2022)

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Golden Gate cloning of epegRNAs

- BsaI-HFv2 (NEB, cat. no. R3733S)
- NcoI-HF (NEB, cat. no. R3193S)
- PvuII-HF (NEB, cat. no. R3151S)
- rCutsmart buffer (10×; NEB, cat. no. B6004S or provided with restriction enzymes)
- Tris-HCl (pH 8.0, 1 M solution; Thermo Fisher Scientific, cat. no. 15568025)
- NaCl (5 M solution; Thermo Fisher Scientific, cat. no. AM9760G)
- T4 DNA ligase (NEB, cat. no. M0202S)
- T4 DNA ligase reaction buffer, 10× provided with the T4 DNA ligase, but can also be ordered separately (NEB, cat. no. B0202S)
- T4 polynucleotide kinase, necessary if sgRNA scaffold oligos for Golden Gate method will be manually phosphorylated (NEB, cat. no. M0201S)
- QIAquick gel extraction kit (Qiagen, cat. no. 28704)

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Component	Amount (μL)	Final concentration
Annealed spacer oligonucleotides, 1 μM Golden Gate part 1 (from Step 3A(iv))	1	0.1 μM
Annealed and phosphorylated sgRNA scaffold oligonucleotides, 1 μM Golden Gate part 2 (from Step 3A(iv) or 3A(vi))	1	0.1 μM
Annealed epegRNA RTT/PBS 3' extension oligonucleotides, 1 μM Golden Gate part 3 (from Step 3A(iv))	1	0.1 μM
pU6-tevopreq1-GG-acceptor 2.2 kb fragment, 30 ng μL <sup>-1</sup> Golden Gate part 4 (from Step 3A(xii))	1	3 ng μL <sup>-1</sup>
BsaI-HFv2 (20 U μL <sup>-1</sup> )	0.2	0.4 U μL <sup>-1</sup>
NcoI-HF (20 U μL <sup>-1</sup> )	0.2	0.4 U μL <sup>-1</sup>
PvuII-HF (20 U μL <sup>-1</sup> )	0.2	0.4 U μL <sup>-1</sup>
T4 DNA ligase (40 U μL <sup>-1</sup> )	0.50	2 U μL <sup>-1</sup>
10× T4 DNA ligase buffer	1	1×
Nuclease-free H <sub>2</sub> O	3.9	—
Total reaction volume	10	—

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TP53-R248Q (c.743G>A)

(参考文献 Nature Biotechnology 2024

High-throughput evaluation of genetic variants with prime editing sensor libraries)

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Spacer:GCAGTCTCTGCCTCGCCAAG (Yellow)

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mExt\_oligo\_top: gtgcTCCAACATGAGTCACTTGGCGAGGC

mExt\_oligo\_bottom: cgcgGCCTCGCCAAGTGA<sup>A</sup>CTCATGTTGGA

NRAS-G12D (c.35G>A) (ClinVar ID: 39648, Accession: VCV000039648.21)

SynDesign (<https://deepcrispr.info/SynDesign>)で設計

Spacer (黄色) Pam (青) PBS-RTT(下線)

Spacer\_oligo\_top: caccGAACTGGTGGTGGTTGGAGgtttt

Spacer\_oligo\_bottom: ctctaaaacCTCCAACCACCACAGTTC

Ext\_oligo\_top: gtgcAACACCATCTGCTCCAACCACC

Ext\_oligo\_bottom: cgcgGGTGGTTGGAGCAGATGGTGTT

Extension Top wSynMut: gtgcAACACCATATGCTCCAACCACC

Extension Bottom wSynMut: cgcgGGTGGTTGGAGCATATGGTGTT

WT sequence:

GTACAAACTGGTGGTTGGAGCAGGTGGTGGTGGGAAAAGCGCACTGACAATCCAGCTAA  
TCCAGAACCACT

Edited Sequence: GGTGGTTGGAGCAG(G/A)TGGTGTT

Edited Sequence wSynMut: GGTGGTTGGAGCA(G/T)(G/A)TGGTGTT