

Construct pLKO.1 shRNA plasmids for gene knockdown

1. Vector backbone

- (1) pLKO.1-TRC cloning vector -Puro (Addgene, #10878)
- (2) pLKO.1-tRFP657 (derived from (1) but exchange Puro into tRFP657)

2. Design shRNA oligos for pLKO.1

Determine 21-mer target in your gene using **Predesigned shRNA** (MERCK)

<https://www.sigmaaldrich.com/JP/en/semi-configurators/shrna?activeLink=productSearch>

-> Search Genes/Clones: key in your interest gene symbol -> search

-> Select shRNA Clones: choose the species and click details. Ex:

<input type="checkbox"/>	TRCN0000148112	Validated	SETDB1	9869	NM_012432	Human	CDS	Details ^
Target Sequence:	GCTCAGATGATAACTTCTGTA			Cell Type:	A549			
Gene Synonym:	ESET, H3-K9-HMTase4, KG1T, KMT1E, TDRD21			Method:	SYBR			
				Knockdown:	95%			
				Version:	1			

Choose validated knockdown efficiency > 90%

Copy Target Sequence as 21-mer as sense, and use complementary sequences as antisense

3. Order Oligos Compatible with pLKO.1

To generate oligos for cloning into pLKO.1, insert your sense and antisense sequences from step 2 into the oligos below. Do not change the ends; these bases are important for cloning the oligos into the pLKO.1 TRC-cloning vector.

Forward oligo:

5' CCGGCC – **21bp sense** – CTCGAG – **21bp antisense** – TTTTGGAAAT 3'

Reverse oligo:

5' AATTATTCCAAAA – **21bp sense** – CTCGAG – **21bp antisense** – GGG 3'

Order oligos with 100 µM DDW

Ex: for human SETDB1

Forward oligo:

CCGGCCC**GCTCAGATGATAACTTCTGTA**CTCGAG**TACAGAAGTTATCATCTGAG**CTTTTTGGAAAT

Reverse oligo:

AATTATTCCAAAA**GCTCAGATGATAACTTCTGTA**CTCGAG**TACAGAAGTTATCATCTGAG**CGGG

4. Cloning Oligos into pLKO.1

The pLKO.1-TRC cloning vector contains a 1.9 kb stuffer that is released upon digestion with EcoRI and AgeI

(1) Digesting pLKO.1 vector

Vector (6 µg)	6 µl
10x rCutSmart	6 µl
EcoRI	1 µl
AgeI	1 µl
DDW	46 µl
<hr/>	
Total	60 µl

Incubate at 37°C for 2 hours and run agarose gel

(2) Annealing oligos

Forward oligo (100 µM)	1 µl
Reverse oligo (100 µM)	1 µl
10x NEB buffer 2	5 µl
DDW	43µl
<hr/>	
Total	50 µl

95°C 4 min

70°C 10 min

Take out the heat block and cool it down to room temperature (at least 1 hour)

(3) Run 0.8% agarose gel (100V, 30 min) with step 4-1 digested vectors

Expected fragment size: 1.9 kb (stuffer), **7 kb (cut vector)**

-> Cut 7 kb digested vector and do gel extraction

-> Elute DNA with 25 µl DDW

(4) Ligation

*10x dilution for annealed oligos: 5 µl mixture (from step 4-2) + 45 µl DDW

10x diluted oligos	3 µl
Purified vector (step 4-3)	1 µl
2x mighty ligation mix	4 µl
<hr/>	
Total	8 µl

Incubate at 16°C for 1-3 hours

(5) Transformation (DH5 α)

- > Take 3 μ l of ligation product (step 4-4) into 25 μ l of DH5 α
- > on ice 30 min -> heat shock 42°C for 45 sec -> on ice 3 min -> add 80 μ l SOC
- > all for plating with Amp plate

(6) Mini-culture

Pick up 4 colonies each group for mini-culture (2 ml LB+2 μ l Amp)

(7) Mini-prep

Dissolve DNA with 30 μ l TE

(8) Enzyme digestion check

Mini-prep DNA	2 μ l
10x buffer M	2 μ l
PvuII	0.3 μ l
DDW	15.7 μ l
<hr/>	
Total	20 μ l

(9) Run 1% agarose gel (100V, 20 min)

Correct expected fragment size: 3712, 2513, 776 bp

Negative control of original vector cutting size: 2810, 2604, 2513, 776, 198 bp

(10) Select correct plasmids for sequencing

*Sequencing method: dGTP

Select Premix1 sequencing and note dGTP 使用 to 特記事項欄

【Premix1 シングル解析】 ¥2,000 / サンプル

テンプレート DNA と Primer を混合したサンプルをご提供いただきます。

シーケンサーにて 2 時間の泳動を行い、ご提供の base 数は 1000bp 程度になります。

Web サイトより Premix1 シングル解析をご選択いただき、「特記事項欄」に、

「dGTP 使用」と予めご記入いただきますと、シーケンス反応時に弊社にて dGTP 試薬を使用して解析を実施します。