


## Viral Infection using patient-derived xenograft (PDX) cells

### ➤ Materials

1. PBS
2. Retronectin (100 ul of 1 mg/ml retronectin diluted with 8 ml PBS)
3. 6-well non-treated plates (or any other small scale)
4. 2% BSA (sterile) in PBS
5. Polybrene (10mg/ml)
6. Virus (Cas9-mCherry as an example)
7. PDX cells

### ➤ Flow chart



➤ **Procedures**

6-well plate as an example (please scale-down if you use smaller plates.)

1. Preparation of diluted retronectin:  
retronectin (1mg/ml, 100 ul/tube) is stored at -30°C common freezer.  
Add 8 ml PBS into a 15 ml tube and then put 100 ul of retronectin into PBS for dilution.
2. Coating retronectin:  
Add 2 ml/well diluted retronectin into a 6-well plate.  
Seal the plate with parafilm to prevent evaporation.  
Put the plate into 4°C fridge for overnight.
3. Coating virus:  
Recycle diluted retronectin into a 15 ml tube (can be reused for 2 times).  
Add 2 ml of 2% BSA in PBS for blocking at room temperature for 30 min.  
Aspirate gently.  
Add 2 ml of PBS for washing.  
Aspirate gently.  
Add 2-4 ml of virus into wells.  
Incubate the plate at a 37°C incubator for 3-6 h.
4. Infection:  
Aspirate virus gently.  
Mix 1 ml of cells (1-2 million) and 1 ml of new virus and 0.8 ul of polybrene stock to make final 4 ug/ml of polybrene.  
\* Add extra cytokines if needed.
5. Spinoculation:  
Put the plate into a centrifuge (attention to balance).  
1000 g, room temperature, 1 h  
Incubate the plate at a 37°C incubator for 24 h.
6. Collect the cells from the 6-well plate with PBS or trypsin:  
Collect suspension cells into a 15 ml tube.  
Add 2 ml PBS into wells and pipet several times to detach the cells from bottom.  
Collect all cells into the same 15 ml tube as the previous step.  
Spin down (1500 rpm, 5 min, 4°C)  
Resuspend the cells with fresh media and perform cell counting.  
Put the cells into an appropriate dish/plate and recover for at least 24 h before further assays.